

## Erlenmeyer Flask Systems

The number one choice of the biopharmaceutical industry in single-use Erlenmeyer Flask assemblies, Bio-Simplex<sup>®</sup> Erlenmeyer Flask Closures are now available solo (Closure Assembly) or with a standard polycarbonate flask (System). Our Bio-Simplex Closure Assembly features the EZ Top<sup>®</sup> assembly, made from pharmaceutical grade C-Flex<sup>®</sup> thermoplastic elastomer (TPE) resin and tubing, and is engineered to ensure an unobstructed fluid path, extremely low levels of extractables and a secure elastomeric seal against plastic, glass and metal surfaces. The Closure Assembly includes C-Flex tubing with a choice of either male Luer connections or CPC<sup>®</sup> quick disconnect fittings, a Bio-Simplex FlexJoint<sup>™</sup>, 0.2 µm vent filter and polypropylene cap. Economically designed for the single use market, both the Bio-Simplex Assembly and System eliminate cleaning and contamination from repeated use. They are easily used with sterile connecting devices and thermal tubing sealers. We offer standard 2 port configurations in 1, 2 and 3 Liter sizes, available as individually packaged Bio-Simplex Closure Assemblies or a complete packaged Bio-Simplex Erlenmeyer Flask System (Bio-Simplex Assembly, Erlenmeyer Flask and cap).

Both come gamma irradiated.


Customer-generated designs, sizes, packaging and other options — such as connectors, heat-sealed ends or alternate materials — are available as custom orders.



*Bio-Simplex Erlenmeyer Flask Systems feature the Closure Assembly, flask and cap.*

EMF SIZE	BIO-SIMPLEX ERLNMEYER FLASK CLOSURE ASSEMBLY (NO FLASK)		BIO-SIMPLEX ERLNMEYER FLASK SYSTEM INCLUDING FLASK
	w/Luer Connection	w/ CPC quick disconnection	w/Luer Connect
1 Liter	EMF1L-LC-EZ-R	EMF1L-QD-EZ-R	EMF1L-LC-R
2 Liter	EMF2L-LC-EZ-R	EMF2L-QD-EZ-R	EMF2L-LC-R
3 Liter	EMF3L-LC-EZ-R	EMF3L-QD-EZ-R	EMF3L-LC-R

## Erlenmeyer Flask Systems



- Closure Assemblies come with either CPC® Quick Disconnect or male Luer connections. Bio-Simplex® Erlenmeyer Flask System comes with standard male Luer end fitting.

- Pharmaceutical grade C-Flex® thermoplastic elastomer (TPE) tubing material is compatible with sterile connecting devices and tube welders.

- Standard polypropylene cap

- Bendable, patent-pending Bio-Simplex Flex-Joint™ retains shape and eliminates kinking.

- Provided double bagged and gamma irradiated.

- 0.2 µm vent filter allows sterile gas exchange.

- Thermal bond between stopper and tubing prevents compromise of product.

- EZ Top® Closure Assembly of high purity pharmaceutical grade C-Flex TPE elastomer provides integral seal to glass, plastic and metal.

- Bio-Simplex Closure Assemblies (minus the flask) are available with either CPC Quick Disconnect or Luer connections in 1, 2 or 3 liter sizes. Bio-Simplex Erlenmeyer Flask Systems (including flask) are available in 1, 2 or 3 Liter sizes.

- Customer-generated designs, sizes, packaging and other options — such as connectors, heat-sealed ends or alternate materials — are available as custom orders.

Bio-Simplex®, C-Flex® and EZ Top® are registered trademarks of Saint-Gobain Corporation. FlexJoint™ is a trademark of Saint-Gobain Corporation. CPC® is a registered trademark of Colder Products Company.

### Saint-Gobain Performance Plastics

4451 110th Avenue North  
Clearwater, FL 33762  
Toll-free: (800) 541-6880  
Tel: (727) 531-4191  
Fax: (727) 530-5603

Leading the way in

**Critical Fluid Transfer and Containment**



**IMPORTANT:** It is the user's responsibility to ensure the suitability and safety of Saint-Gobain Performance Plastics products for all intended uses and that the materials to be used comply with all applicable medical regulatory requirements. Saint-Gobain Performance Plastics assumes no responsibility for any product failures that occur due to misuse of the materials it provides arising out of the design, fabrication or application of the products into which the materials are incorporated.

SAINT-GOBAIN PERFORMANCE PLASTICS DISCLAIMS ALL IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

WARRANTY: For a period of 12 months from the date of first sale, Saint-Gobain Performance Plastics warrants this product to be free of defects in materials and workmanship. Our only obligation will be to replace any portion proving defective, or at our option, to refund the purchase price thereof.